

Sample Printout  
using  
\*NHANES III  
predicted values.


\*NHANES III predicted values  
come standard on the Alpha III.  
Other values available.

VITALOGRAPH - ALPHA

DATE: 10/ 2/03  
 NAME: .....  
 REFERENCE NO: 1234  
 AGE: 61  
 HEIGHT: 60 INCHES  
 SEX: FEMALE  
 ETHNIC ORIGIN: CAUCASIAN  
 PREDICTION REFERENCE: HANKINSON  
  
 SOFTWARE REF.NO: 61257 ISSUE 1  
 CALIBRATION UPDATE DATE: 8/22/03  
 CALIBRATION CHECK DATE : 10/ 2/03

PRE : 3 TESTS, BEST WITHIN 3%  
 POST: 3 TESTS, BEST WITHIN 2%

	PRE			POST			CHANGE
	PRED	MEAS	%	MEAS	%	%	
VC	2.80	----	-	----	-	-	-
FVC	2.80	2.59	92	2.59	92	0	0
FEV1	2.16	1.85	86	1.87	87	1	1
FEV1 %	78	71	-7	72	-6	1	1
PEF	339	306	90	300	89	-1	-1
FEF25-75	2.08	1.09	52	1.20	58	10	10
FIF25	3.53	2.50	71	2.60	74	4	4
FIF50	3.82	2.50	65	2.50	65	-	-
FIF75	2.60	2.80	108	2.40	92	-15	-15



**\*\* INTERPRETATION \*\***  
 NORMAL VENTILATORY FUNCTION  
 SLIGHT IMPROVEMENT POST-BRONCHODILATOR  
 CURVE VERIFICATION: SIGNED .....