



the difference is clear

Vitalograph BreathCO

Smoking Cessation Motivation

- ❖ Measures breath CO ppm
- ❖ Hand-held, battery operated
- ❖ Quick and easy to use
- ❖ Non-invasive
- ❖ Easy to read, large display
- ❖ Unique, accurate zero adjuster slide
- ❖ Reliable calibration
- ❖ Inexpensive, hygienic operation
- ❖ Excellent for individual or group therapy

The heart of successful smoking cessation is motivation. Improved health is an excellent motivator in getting people to kick the smoking habit for good.

The Vitalograph BreathCO provides an easy, noninvasive, and immediate way of assessing a participant's smoking status and provides visual proof of their breath CO concentration. Teaching them how this reflects the dangerous levels of poisonous CO in their lungs that enters into their blood stream and carried throughout their bodies, lowering the blood oxygen, shows them the damage they are doing to themselves as a result of the habit, especially the heart.

Incorporating the BreathCO as a motivator as well as an indicator of smoking status with your arsenal of other smoking cessation tools and therapy creates a well balanced smoking cessation program that is bound for success.





the difference is clear



Technical Specification:

Product:

BreathCO monitor

Model Number:

29700

Size:

7" x 4" x 2"

Weight:

10 oz

Operating Temperature:

32 - 100 °F

Storage Temperature:

25 - 110 °F

Measuring Accuracy:

± 3 ppm

Display Resolution:

1 ppm

Display Range:

0 to 199 ppm

Linear Range:

0 to 99 ppm

Humidity Range:

0 to 99%

Power Supply:

9V Dry Alkaline Battery

Digital Display:

3.5 digit 25mm LCD

Ordering Info:

Catalogue No.

29700 BreathCO Monitor

29040 Mini mouthpieces (40)

20202 Mini mouthpieces (200)

29512 T-Valves (10)

29504 CO Calibration Kit

29506 CO Calibration Gas Refill

Returning the BreathCO annually for recalibration and service is recommended.

Vitalograph Inc.
13310 West 99th Street
Lenexa
Kansas 66215 USA
Tel.: (913) 888 4221
Fax: (913) 888 4259
E-mail: vitcs@vitalograph.com

www.vitalograph.com

